ELITE 900

ENGLISH

Upper Arm Blood Pressure Monitor

Instructions for use3

ESPAÑOL

Tensiómetro de brazo

Manual de instrucciones..27

FRANÇAIS

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Mode d'emploi47

For Canada only/Solo para Canadá/
Pour le Canada seulement :
Manufacturer/ Fabricante/ Fabricant:
Beurer GmbH, Söflinger Str. 218
89077 Ulm Germany

Distributed by/Distribuido por/
Distribué par :
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.shop-beurer.com
Assembled in China

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**READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE
USING THIS PRODUCT**

Keep this manual in a safe location for future reference

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ENGLISH

Contents

1. Important safety notes	4	8. Displaying and Deleting Measurements	14
2. Intended Use	8	9. Care, Maintenance and Disposal	15
3. Package Contents	8	10. Troubleshooting Guide	16
4. Parts and Controls	8	11. Technical Specifications	17
5. Preparing for Measurement	9	12. FCC Compliance Information	18
6. Measuring Blood Pressure	11	13. Warranty	19
7. Evaluating Results	13		

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact customer service. See warranty for service contact.

Dear Valued Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for application in the areas of heating, weight management, blood pressure, body temperature, pulse, gentle therapy, massage, and beauty.

Register your product today to receive further benefits:

- added protection for your purchase
- health and wellness tips
- new product information
- continual outstanding customer service and support

Visit www.registerbeurer.com or contact customer service at **1-800-536-0366**.

Please read these instructions for use carefully and keep them for later use.

Make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer team

1. Important safety notes

Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
WARNING	Indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.
CAUTION	Indicates a potentially hazardous situation, which, if not avoided, may result in minor/moderate injury.
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage.
	Product information Note on important information
	Follow instructions for use

	Type BF, applied part
	Direct current
Storage/Transport 	Permissible storage temperature and humidity
Operating 	Permissible operating temperature and humidity
	Atmospheric pressure, limit Indicates the range of atmospheric pressures to which the medical device can be safely exposed
IP22	Device protected against foreign objects ≥ 12.5 mm and against water dripping at an angle
	Serial number
	Protection class II device The device is double-insulated and, therefore, is in protection class 2

WARNING

- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. Consult your doctor before using the blood pressure monitor during pregnancy.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who has responsibility for their safety or they receive instructions from this person on how to use the device. Supervise children around the device to ensure they do not play with it.

- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders, and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Using the blood pressure monitor outside your home environment or whilst on the move (e.g. whilst travelling in a car, ambulance or helicopter, or whilst undertaking physical activity such as playing sport) can influence the measurement accuracy and cause incorrect measurements.
- People with pacemakers or other electrical implants should consult their doctor before using the device.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.
- Do not perform measurements more frequently than necessary. The resulting restriction of the blood flow may cause injury.
- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.
- Small parts may present a choking hazard for small children if swallowed. They should therefore always be supervised.

CAUTION

- Remove any kind of upper arm jewelry or the like before taking a measurement. This could cause bruises.
- Do not use the device at the same time as other medical electrical devices (ME equipment). This could lead to a malfunction of the device and/or an inaccurate measurement.
- Do not use the device if you have metal implants.

NOTICE:

- Never immerse and or spill water or any other liquid onto the monitor or any components, otherwise liquid will enter it and cause damage.
- Never attempt to repair, open and/or disassemble the unit (including arm cuff and optional accessory) or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact.
- Bring the device to room temperature before taking measurements. If the measuring device has been stored near to the maximum or minimum storage and transport temperature and is brought into an environment with a temperature of 68 °F (20 °C), it is recommended that you wait 2 hours before using the measuring device.
- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
 - Only use the cuff included with the delivery or original replacement cuffs. Otherwise incorrect measured values will be recorded.

Battery Handling Safety Precautions

- Do not remove or attempt to remove the non-user-replaceable battery.
- Risk of explosion! Never throw batteries into a fire.
- Batteries must not be short-circuited.
- Only use chargers specified in the instructions for use.
- Fully charge the battery prior to initial use and at least every 3 months.
- Protect batteries from excessive heat.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.
- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- Do not disassemble, open or crush the batteries.

2. Intended Use

This blood pressure monitor performs non-invasive measurement and monitoring of arterial blood pressure and heart rate on adult patients.

3. Package Contents


1 x Upper Arm Blood Pressure Monitor with Cuff (Cuff size: 8.7 - 16.5 in / 22 - 42 cm)

1 x USB-C Cable

1 x Quick Start Guide


1 x Instruction manual

4. Parts and Controls

- 1 LED display
- 2 **MEMORY** Button 
- 3 **START/STOP** Button 
- 4 USB-C Charging Port
- 5 Arm Cuff

LED Display Information:


- 6 Bluetooth® Indicator ✱
- 7 Systolic Pressure
- 8 Unit of Measure
- 9 Diastolic Pressure
- 10 Pulse Icon ♥
- 11 Display of Values
- 12 Risk Indicator Scale
- 13 Irregular Heartbeat Icon

5. Preparing for Measurement

5.1 Initial use

Charge the batteries:

- Charge the device fully before using.
- If  appears on the display, you must charge the device.
- Connect the USB-C charging cable (included) to the USB-C charging port and to a USB power supply or computer.

- Press the **START/STOP** Button (ⓘ) or **MEMORY** Button (📖) at any time, 🔌 blinks while the battery is charging. 🔌 appears once the device is fully charged.
- The display switches off after a few seconds.

The device cannot be used while it is being charged.

Before using the device for the first time, you will need to pair it with the “beurer HealthManager Pro” app in order to synchronize the time and date. To do so, proceed as follows:

5.2 Bluetooth® Setup

- Download the free “beurer HealthManager Pro” app from the Apple App Store or Google Play.

Click here for the
“beurer HealthManager Pro” app

List of system requirements
and compatible devices

- Activate *Bluetooth®* in your smartphone’s settings.
- Start the app.
- Select ELITE 900 in the app and follow the instructions.

5.3 Connect to the “beurer HealthManager Pro” app

For the device to store the correct time and date, it must connect to the mobile app periodically. Open the “beurer HealthManager Pro” app on your mobile device and turn on the blood pressure monitor by pressing the **START/STOP** Button (ⓘ).

- If the device detects that the time has not been set, “APP” will show on the LED Display.
- No successful *Bluetooth®* connection within 30 seconds: “NO TIME” is shown on the LED display.
- *Bluetooth®* connection successful: Time is being synchronized
- After the time is synchronized: “TIME” appears on the LED Display.

To skip setting the time and date, first press the **START/STOP** Button (①) to turn on the device. “APP” will show on the display. Press the **START/STOP** Button (①) again and “NO TIME” will show on the LED display. Note that any subsequent measurement will have an incorrect time and date associated with it in the mobile app.




6. Measuring Blood Pressure

General rules when measuring your own blood pressure:

- The blood pressure monitor is suitable for home use and needs to be at room temperature when measuring.
- Measurements can be performed on either arm. For consistency, always measure your blood pressure on the same arm.
- For consistency, always measure your blood pressure at the same time each day.
- When performing multiple measurements on an individual, wait at least 5 minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Before the initial blood pressure measurement, make sure always to rest for about 5 minutes.
- Repeat the measurement if you are unsure of the results.



Positioning the Cuff:

1. This unit includes a universal size cuff that fits upper arms 8.7 in – 16.5 in (22 – 42 cm) in circumference. Place the cuff round your bare upper arm 0.8 - 1.2 inches (2-3 cm) above the bend of the elbow and above the artery. The center of the device should form a line with the center of the palm (see picture). Ensure blood circulation is not restricted by tight clothing or other objects.
2. Fasten the cuff, making sure it is not too tight around your arm.
3. You should be able to fit two fingers under the cuff after it is in place.


Correct Posture:

- Rest for at least five minutes before each measurement.
- Perform the measurement either sitting or lying down. Always make sure that the Cuff is level with your heart.
- To perform a measurement in a seated position, sit comfortably with your arms and back supported. Do not cross your legs and place your feet flat on the floor.
- Keep still during the measurement and do not talk.
- Relax your arm and the palm of your hand.

Performing a Measurement:

- With the cuff attached properly and in a correct posture, press the **START/STOP** Button (①). All LED Display elements will briefly display and measurement will begin. ♥ is displayed as soon as a pulse is detected. To stop a measurement at any time, press the **START/STOP** Button (①).
- As a measurement concludes, air pressure in the cuff will release. Systolic/diastolic pressure (along with risk indicator) and pulse rate will display alternately.

- “Er_” (where the blank space is a number 1-7) will appear if the measurement was not performed properly. Refer to the “Troubleshooting Guide” section for details about specific error message numbers and then repeat the measurement.

The device switches off automatically after 30 seconds if no button has been pressed or press the **START/STOP** Button (①) consecutively for 2 seconds to switch off the device.

Transfer of measured values via *Bluetooth*®

- Make sure the “beurer HealthManager Pro” app on your smartphone is activated and open before turning on the blood pressure monitor.
- To transfer the measured values via *Bluetooth*®, the device first connects to the “beurer HealthManager Pro” app. The *Bluetooth*® icon (✱) will flash during this process.
- Once the device has successfully connected to the app, the *Bluetooth*® icon (✱) will be shown permanently on the display.
- If the device cannot establish a *Bluetooth*® connection, the device will shut down automatically after about 30 seconds.
- With a successful *Bluetooth*® connection, measured values are transferred automatically.
- If data transfer is unsuccessful, the error message “Er7” is displayed.

7. Evaluating Results

Irregular Heartbeat:

This device can detect irregular heartbeat patterns during measurement and will indicate them by displaying the Irregular Heartbeat Icon . This could be an indicator for arrhythmia, a condition where the heart rhythm is abnormal. These symptoms (omitted or premature heartbeats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress, or lack of sleep. Only your doctor can determine if you actually have arrhythmia.

Repeat the measurement if the Irregular Heartbeat Icon displays after a measurement. Rest for at least five minutes between measurements and do not talk or move during measurement. If the icon appears often, contact your doctor, since self-diagnosis and treatment based solely on the test results may be dangerous. It is vital to follow your doctor's instructions.

LED Risk Classification

Measurements can be classified and assessed by the table below. However, these standard values are only a general guideline, since individual blood pressure varies in different people and different age groups, etc.

Category	Systolic (in mmHg)	Diastolic (in mmHg)	Risk Indicator
Grade 3: Severe hypertension ¹	≥ 180	≥ 110	Red
Grade 2: Moderate hypertension ¹	160–179	100–109	Orange
Grade 1: Mild hypertension ¹	140–159	90–99	Yellow
High Normal ¹	130–139	85–89	Green
Normal ¹	120–129	80–84	Green
Optimal ¹	< 120	< 80	Green
Low blood pressure, hypotension ²	< 90	< 60	Orange

¹Source: WHO, 1999 (World Health Organization)

²Source: American Heart Association, 2023

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous. The classification on the display shows which category the recorded blood pressure values fall into. If the systolic and diastolic values fall into two different categories, the Risk Classification Indicator on the device always shows the higher category.

Low blood pressure

A blood pressure measurement of less than 90/60 mmHg is considered hypotension (blood pressure that is too low). A single lower-than-normal reading is usually not cause for alarm unless you are experiencing any other symptoms or problems. Seek medical attention if you experience low blood pressure accompanied by dizziness, lightheadedness, nausea or other symptoms. Self-diagnosis and treatment based solely on measurement results may be dangerous.

8. Displaying and Deleting Measurements

Every successful measurement is stored along with its date and time. If there are more than 240 measurements in the memory bank, the oldest measurement is deleted.

Displaying Measured Values:

With the device shut off, press the **MEMORY** Button () to display the stored measurements in order, starting with the latest measurement.

Press the **MEMORY** Button () again to display the previous measurement or press the **START/STOP** Button () to display the next measurement.

To turn the device off, press and hold the **START/STOP** Button () for about two seconds.

NOTE: If the time and date have not been set via the “beurer Health Manager Pro” app, “APP” will be shown on the LED display after pressing the **MEMORY** Button (). To skip setting the time and date via the app, press the **START/STOP** Button () and “NO TIME” will be displayed. Press the **MEMORY** Button () again to view the stored measurement values as described above.

Clearing Memory:

Access the measured values as explained above, and press and hold the **MEMORY** Button () for two seconds; “no” will appear in the LED Display and all values in memory will be deleted.

9. Care, Maintenance and Disposal

- Clean the device and cuff with a slightly moistened cloth.
- Do not use detergents or solvents.
- Never immerse the device in water or other liquids.
- When storing the device, make sure no heavy objects are placed on top of it.

Disposal

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

10. Troubleshooting Guide

Error message	Possible cause	Solution
ER1	Pulse not detected.	Please wait one minute and repeat the measurement. Ensure that you do not speak or move during the measurement.
ER2	Measured values are outside the measurement range.	
ER3	Cuff is attached incorrectly.	Refer to the “Positioning the Cuff” section for details.
ER4	An unknown error occurred during measurement	Please wait one minute and repeat the measurement. Ensure that you do not speak or move during the measurement. If the error occurs repeatedly, consult a doctor to check you are healthy.
ER5	Cuff does not inflate properly. (pressure is higher than 300 mmHg)	Please take another measurement to check whether the cuff can be correctly inflated.
ER6	System error.	Contact Customer Service.
ER7	There is a problem with the connection between the device and the app.	Turn off the device, close the app, and deactivate Bluetooth on your mobile device. Reactivate Bluetooth and try to establish the connection again.
	The battery is almost empty.	Charge the battery before using again.

11. Technical Specifications

Model	ELITE 900
Type	BM 59
Measurement Method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement Range	Cuff pressure 0–295 mm Hg, systolic 57 – 255 mm Hg, diastolic 25–195 mm Hg, Pulse 40–199 beats per minute
Display Accuracy	Pressure ± 3 mmHg, Pulse ± 5 %
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory	1 x 240 memory spaces
Dimensions	4.9 in L x 1.9 in W x 1.1 in H (125 mm L x 48 mm W x 28 mm H)
Weight	7.9 oz (225 g) without batteries
Cuff size	8.7 in - 16.5 in (22 - 42 cm)
Operating Range	41° F - 104° F (5° C - 40° C), 15-90% relative humidity (non-condensing) Atmospheric pressure range of 700 hPa to 1060 hPa
Storage /Transport Range	-4° F - 131° F (-20° C - 55° C), 10 - 90% relative humidity (non-condensing)
Power	Charge: 5V $\overline{\overline{\overline{\quad}}}$ 1A Li-ion rechargeable battery, 3.7 V

Number of measurement cycles
per battery charge

Approx. 60 measurements

The serial number is located on the device.

This device contains licence-exempt transmitter(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

(1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device. This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

12. FCC Compliance Information

Upper Arm Blood Pressure Monitor ELITE 900

Responsible Party – U.S. Contact Information

Beurer North America LP, 1 Oakwood Boulevard, Suite 255, Hollywood, FL 33020. United States

1-800-536-0366, info@beurer.com

FCC Compliance Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio

communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

13. Warranty

Your Beurer Upper Arm Blood Pressure Monitor, model ELITE 900, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Upper Arm Blood Pressure Monitor, model ELITE 900, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage;

theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you. Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

For more information regarding our product line in the USA, please visit: www.shop-beurer.com

For Canada only:

Manufacturer:

Beurer GmbH

Soeflinger Str. 218

89077 Ulm, Germany

Distributed by:

Beurer North America LP

1 Oakwood Boulevard, Suite 255

Hollywood, FL 33020, USA

www.shop-beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

Assembled in China

EMC Guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home environments.

Warning: Don't use near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment , including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1**Guidance and manufacturer's declaration - electromagnetic emissions**

Emissions test	Compliance
RF emissions, CISPR 11	Group 1,
RF emissions, CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions, IEC 61000-3-3	Not applicable

Table 2**Guidance and manufacturer's declaration - electromagnetic Immunity**

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable

Immunity Test	IEC 60601-1-2 Test level	Compliance level
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	Not applicable	Not applicable
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 3**Guidance and manufacturer's declaration - electromagnetic Immunity**

	Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27	27
	450	430-470	GMRS 460, FRS 460	FM $\pm$ 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9	9
	745							
	780							
	810	900-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28
	870							
	930							

	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
	1845							
	1970							
	2450	2400-2570	Bluetooth®, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9
	5500							
	5785							

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The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.

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