

4-WHEEL AIR FRAME WALKER

MOB1117



To see all the FAQ's,
visit [vhealth.link/t3z](https://vivehealth.link/t3z)

vive
vivehealth.com

QUICK START GUIDE

FOLDING AND UNFOLDING THE WALKER:

To unfold the Walker, pull apart the frame by the two handles. The lock will click into place once fully unfolded.

Locked Walker:

To fold the Walker, Pull up on the locking lever on the right side of the frame. Then continue to push the frame together.



Moving the Handles up and down:

Unscrew the star knob screw on the front uprights. Once removed you may set the Handles to the desired height and then re-insert the star knob screw.



USING THE BRAKES:

Pull up on the brake handle to engage the brakes.

Note: This also locks the brakes.



vive
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

Releasing the brake lock:

Pull up on the smaller red lever underneath the brake handles.



WARNINGS ⚠

- 300lbs weight capacity
- Do not use strong corrosive cleaner to clean the product's surface.
- Product should not be left outside in rainy conditions. Product may rust if regularly exposed to moisture.
- Do not hang any heavy weighted items on either of the armrests as it will cause the Walker to be unstable.

- This item is meant to be a walking aid. This means that it is not meant to support the full weight of one's body, it is only meant to assist you while walking.
- Make sure the device is fully opened and clicked into the "locked" position before use.
- Regularly check the device for any broken or cracked pieces or parts.
- Vive is not responsible for injury or death caused by improper or unintended use.