

MATERNITY BELT

SUP3085



To see all FAQs in one place,
visit [vhealth.link/rcm](https://vivehealth.link/rcm)

vive
vivehealth.com

QUICK START GUIDE

HOW TO USE

Note: Depending on the size of your baby bump, assistance may be needed to put this on.

1. **Position the Abdominal Support Pad:**

Place the pad under your belly and secure it by wrapping the strap around your body. Make sure it's low as possible without blocking your legs when you walk or sit.



2. **Attach the Back Support Brace:**

Connect one end of the brace to the abdominal pad, wrap it around your back, and secure it to the other side.



3. **Secure the Elastic Tummy Strap:**

Attach one end of the strap to the back brace. Stretch it lightly over your belly and secure it to the other side to keep the belt in place.



vive
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

4. **Adjust the Back Support:**
Pull the elastic straps on each side of the back brace to customize the support and comfort.



CARE INSTRUCTIONS

- Hand wash in cold water with mild soap.
- Air dry flat.
- Do not bleach, iron, or dry clean.
- Store in a cool, dry place.

PROPER FIT & USAGE

- Belt should feel snug but not tight.
- Support pad sits under belly.
- Wear up to 12 hours daily.
- Remove if you feel:
 - Discomfort
 - Pinching
 - Restricted breathing

WARNINGS ⚠

- Consult your healthcare provider before use.
- Discontinue use if discomfort occurs.
- Do not wear while sleeping.
- Not intended to diagnose, treat, cure, or prevent any disease.
- Consult your healthcare provider if you notice any changes in your condition.