

CORE SHOWER COMMODE TRANSPORT CHAIR

MOB1115

vive®
MOBILITY



To see all FAQs in
one place, visit
visit vhealth.link/h3y

vivehealth.com

CORE SHOWER COMMODE TRANSPORT CHAIR

WHAT'S INCLUDED



- Frame (A)
- 2x Rear Legs (B)
- 2x Front Legs (C)
- 2x C-clips (D)
- Seat (E)
- Bucket (F)
- Backrest (G)
- Footrest (H)

INSTRUCTIONS



For a video demonstration,
check out vhealth.link/m4m

1. Unfold the Frame portion of the Commode Chair. Snap the front bar portion into the openings on each side of the vertical bars that swing around to the front. Make sure they click into place with the locking pins.
2. Add the 2 Rear Legs with Wheels (B) to the bottom of the frame. The rear legs will match up with the posts on the frame that have a gold pin sticking out for the rear legs to click into, note there are 3 height settings make sure all 4 legs are installed at the same height.



3. Install the front wheels (C), these will be secured with the (D) C-clip pins. Again, Install these at the same height setting as the back legs
4. Install the Seat and the hole cover for the commode portion to the frame. The seat will clip onto the rear bar and rest safely on the front bar.
5. On the Frame (A) Attach the backrest so the pins click into the provided holes of the backside of the frame.
6. Insert the commode bucket (F) into the rear slot under the seat.



7. Insert the Footrest (H) onto the two pipes that jet out from the two front legs.



Use the C-Clips (D) to secure the Footrest in place.



8. Underneath the footrest, There is flip down legs that work as stoppers. These legs are for use during prolonged duration in the chair like showering. They work excellent to help hold the chair in place especially on slippery surfaces.



LOCKING THE BRAKES

Push down on the locking tab with your foot on all 4 wheels to properly and safely lock the brakes. Push upward with your foot on the same tab to unlock them.



CARE INSTRUCTIONS

- The chair can be wiped down with mild soap, an ordinary bathroom surface cleaner, or disinfectant.
- Best to wipe dry after use in the shower.

WARNINGS ⚠

- Maximum Weight Capacity: 275LBS
- Never try to stand up from the chair on a wet surface.
- Never try to stand up from the chair while the brakes are not engaged.
- This wheelchair has not been approved as a seat that can be used within vehicles of any kind.

- Do not attempt to operate chair on inclines without proper assistance.
- Do not operate this wheelchair on roadways, streets or any surfaces with vehicular traffic.
- Do not stand on a wheelchair.
- Abruptly changing directions while going down an incline could cause instability and/or increase the possibility of tipping over.
- Wheelchair should be on a smooth, stable and level area with the hand brakes engaged by the transporter before transferring the passenger in or out of the wheelchair. When transferring in or out of a wheelchair, do not stand on or apply weight to foot rests.
- Do not lean over or apply weight to the top of the wheelchair back. This can cause the wheelchair to tip over.

CAUTION:

Please note that this chair is designed to be used for transporting another individual and is not meant for self transportation. This Transport Chair is designed to safely be a stable center of gravity but it is important to try and keep a stable center of gravity at all time while performing activities that may change your center of gravity or a passenger's center of gravity. Never reach out so far that it would require you to shift in the seat or bend over out of the chair. When you have to reach for objects that require you to lean, be certain your front wheels are pointing forward as this will stabilize the chair.

Ramps: Be sure of your own capabilities, limitations and strength before attempting to push another individual up a ramp. Correct techniques on ramp negotiation should be shown to you by a qualified professional. If you must stop on an incline, avoid any and all sudden or abrupt movement when you resume movement. Because of the backward angle of the chair, sudden abrupt movement could cause instability. Do not attempt to change direction when moving in incline or decline, this could result in tipping.



HAVE MORE QUESTIONS?

Check out our list of Frequently Asked Questions at vhealth.link/h3y for helpful answers.



And if that doesn't answer your question, our customer service team would love to help! Feel free to connect with them by phone, e-mail, or chat on our website

service@vivehealth.com
1-800-487-3808
vivehealth.com

Distributed by



8955 Fontana Del Sol Way
Naples, FL 34109
Made in China