

Thera **Lite**[®]

Aura

Mood and Energy Enhancing Light

USER GUIDE

Thank you for purchasing the TheraLite Aura.

Please read this manual in its entirety. For more information on light therapy, please visit carex.com.

Theralite Aura has been designed with your needs in mind. Projecting light from above, adjustable angle and multiple light settings make the TheraLite Aura the perfect companion at home or at the office.

UV-free LED technology and a comfortable color temperature ensure you experience a safe, comfortable therapeutic light. The adjustable lamphead lets you switch from therapy lamp to desk light with the turn of the knob.

Bright Light Benefits

Daily exposure to bright light helps to maintain your circadian rhythm, boost energy and manage sleeping patterns.

For many, morning sun exposure is limited by work, school or other commitments, even in the summer. The TheraLite Aura is a bright light therapy lamp designed to fill this gap.

INCLUDED ITEMS

TheraLite Aura and User Guide



⚠ WARNING: This light should not be used on patients with macular degeneration. **ALWAYS** consult a physician before using this product.

TheraLite Aura Light Therapy Desk Lamp
Model Number: FGP83000 0000

OPERATING INSTRUCTIONS

⚠ WARNING: RISK OF ELECTRICAL SHOCK

- Turn off the power before performing any maintenance.
- Verify that supply voltage is correct by comparing it with the lamp label information.
- If you have eye problems or have had eye surgery, you should consult your doctor or ophthalmologist before you use any form of light therapy. Certain eye and eye-related conditions may respond unfavorably to light therapy.
- As with any electrical device, **DO NOT** use the TheraLite in or near water.

⚠ CAUTION: RISK OF INJURY

- Avoid direct eye exposure to the light source while it is on.

General Safety Information and Precautions

- This product is suitable for indoor use **ONLY**.
- **DO NOT** overuse the TheraLite. Overuse may cause irritability, excessive energy or difficulty falling asleep at bedtime.
- It is **NOT** necessary to stare into the light. Read, work or talk on the phone during your light therapy session.
- During the first few uses, your TheraLite may give off an odor; this is harmless, will **NOT** affect the use, and will quickly disappear.
- **ALWAYS** use in a well-ventilated area away from other heat sources and ensure that air can flow freely around the lamphead when the product is in use.
- To reduce the risk of death, personal injury or property damage from fire, electric shock, falling parts, cuts/abrasions and other hazards, read all warnings and instructions included with and on the fixture box and all fixture labels.
- For residential installation: If you are unsure about the installation or maintenance of the lamps, consult a qualified licensed electrician and check your local electrical code.
- **DO NOT INSTALL DAMAGED PRODUCT!**
- The external flexible cable or cord of this lamp cannot be replaced; if the cable or cord is damaged, the lamp will be destroyed.
- To avoid injury or damage to the fitting, ensure that power leads and screws are secure before connecting the power.
- Young children should be supervised to ensure that they **DO NOT** play with the portable lamp.
- Handle the fitting and/or shade with care to avoid any damage or breakage.
- Select a suitable location away from liquids and hazards.
- Ensure that the fitting **DOES NOT** come in contact with corrosive chemicals, etc.
- To clean, wipe with a damp clean cloth. **NEVER** soak the fixture with water.

- ⚠ WARNING:** This product emits a bright light which may contribute to migraines in some individuals.

Installation Requirements


1. Select a suitable location for the lamp.
 - Ensure it is on a flat and stable surface.
 - This product is suitable for indoor use **ONLY**.
 - **DO NOT** place flammable material near or over the fixture.
 - Ensure the fixture is not installed next to any air vents and is **NOT** in any moist or damp areas.
2. Take care **NOT** to pull any electrical wires during unpacking as this may damage the connection.
3. Check whether the fitting has been damaged during transport.

DO NOT operate/install any product which appears damaged in any way. Return the complete product to the place of purchase for inspection, repair or replacement.

DO NOT plug into the main power socket outlet before or during assembly. Only do so after the portable lamp is fully assembled and ready to use for the first time.

Installation Directions

1. Remove all packaging material from the product.
2. Select a suitable location for the light fitting. See *Installation Requirements* above.
3. Adjust the head of the lamp to the desired angle.
4. Insert the fixture plug into a properly grounded outlet. Ensure it is fully plugged in.
5. Use switch to turn the lamp on/off.




Assembly Instructions

1. On a level, solid surface lay out the lamphead/arm, the lamp base and the 2 Philips screws for assembly.
2. Position the lamphead/arm on its side.
3. Carefully and firmly insert the bottom of the lamphead/arm into recess on the lamp base in the orientation shown.
4. Insert and by hand seat the Philips screws through the holes on the bottom of the lamp base.
5. Using a Philips screwdriver, hand-tighten the screws until snug. Be careful not to overtighten.
6. Stand the lamp upright.
7. While holding the lamphead with one hand, with the other, turn the angle adjustment knob to loosen the lamphead, position to the desired angle, and re-tighten the adjustment knob.
8. Plug the lamp into a standard 120V electrical outlet.

Operating the Light

1. Use the Light ON/OFF [⏻] touch button to turn the lamp on.
2. Press ⊕ to increase the brightness. Press ⊖ to lower the brightness.



Using Your TheraLite for Therapy

For those with sleep or circadian rhythm issues, regular daily use of the TheraLite Aura can help provide relief. For most conditions, use shortly after waking or in the morning.


Place the TheraLite Aura above your line of sight, where possible, to maximize the light reaching your eyes. You should be facing the lamp with eyes open during sessions, but you do not need to look directly at the light. Over the first few days of use, assess the effect and, if necessary, lengthen or shorten the session time based on your response.

Symptoms of Overuse

Negative effects are rare with light therapy. However, if you experience increased irritability, excessive energy, and/or any consistent discomfort during or after use, decrease session time or move the light further away. Using the lamp late in the afternoon or in the evening is not advised, as this may disrupt your sleep. It is recommended that you consult your doctor before beginning any bright light regimen.

Light Therapy Basics

LUX or the light intensity at a specific distance from the light is the common term used in bright light therapy. The chart on the right details the average treatment times based on your eyes' proximity to the lamp.



For Therapy: Adjust to the best angle so the light shines downward toward your eyes.

DISTANCE	SESSION TIME
12" / 30cm	30 minutes
18" / 45cm	60 minutes
24" / 60cm	120 minutes

Specifications

Power: 120V, 60Hz, 35W

Size: 22" x 11.25" x 11.5" / 56cm x 28cm x 29cm

Lens: High-impact polycarbonate

Light Intensity: 10,000 LUX at 12" / 30cm

Light Source: 168 LEDs, 3000K, 100% UV Free

Angles: Adjustable

Warranty

Your lamp is maintenance-free and should provide years of service. The TheraLite Aura has a one-year limited warranty from the date of purchase by the original purchaser. Please retain proof of purchase to ensure eligibility for warranty coverage. Compass Health Brands may choose to repair or replace a defective unit. Shipping costs are the responsibility of the product owner. Misuse, abuse, accidental or incorrect operation, contrary to the instructions, will void a warranty claim.

COMPASSHEALTH

Manufactured for: Compass Health Brands Corp.
6753 Engle Rd., Middleburg Heights, OH 44130 USA
www.compasshealthbrands.com
800.376.7263 • Made in China
42-FGP83000 0000_04 © 2025



Disclaimer

TheraLite is an innovative light supply system and is not a listed medical device in the USA.

READ USER GUIDE PRIOR TO USE.