


A fresh flossing adventure



Welcome, bacteria slayer

to our elite order of fearless flossers! This kit will give you everything you need. Your quest? Protect those pearly whites and vanquish plaque wherever it may hide.

The fate of the realm (ahem, your mouth) rests in your hands.

“Go forth and get flossing!”

Hamish 😊 Brittany 😊

Hamish & Brittany
Best friends and founders of BURST

Prep your case

1. To open, turn over the floss case and squeeze the two hooks on the underside.
2. Pull off the floss cover.
3. Place the bobbin on the spool and thread the floss around the dispensing post.
4. Replace the floss cover and pull a good amount through the cutter.
5. Close the lid.



How to floss

By Charlotte Stirk, Registered Dental Hygienist

Step 1

The setup

Gather a generous length of floss around one middle finger and wind a few inches around the other.

Step 2

Get a grip

Pinch the floss with your thumb and pointer fingers.

Step 3

Glide


Guide the floss between every tooth, shimmying up and down toward the gumline. Make a C-shape so it can get right up under those gums!

Step 4

Keep it moving

As you clean each chomper, wind used floss around the middle finger you've designated the “dirty” section so the entire mouth gets the fresh floss treatment.

Check out our video tutorial from a real dental hygienist!



Floss refills

Already subscribed? You're all set! Expect a fresh bobbin in the mail every 12 weeks. If you haven't subscribed yet, click below to get started and keep that smile on schedule.

[Subscribe to floss now](#)

The dream clean team

Just like every BURST product, these mouth-saving superstars were made with our hygienist and dentist Ambassadors.



Toothpaste

Looking to crush plaque, fight tooth decay, and freshen breath? You've come to the right paste.

[Buy now](#)

Sonic Toothbrush

With three marvelous modes, this brush gives teeth and gums the royal treatment – smiles never looked (or felt) so good!

[Buy now](#)



Water Flosser

Spray away plaque, freshen breath, and massage gums for a powerful, elevated pampering experience.

[Buy now](#)



Questions? We've got the answers

Is the BURST Floss waxed?

Yes! The BURST Floss uses Microcrystalline wax, there are No PFAs.

Why does the BURST floss expand and how?

The BURST Floss is made with technology that reacts to the saliva and friction in your mouth, which then causes the filament to expand. The expansion of the filament will reach into the nooks and crannies in your teeth to get rid of your plaque.

Can I use the BURST Floss with implants?

We would suggest that you consult with your doctor regarding the appropriate at-home care for your implants.

How is the BURST Floss different than other Floss?

Our BURST Floss is beautifully packaged in an eco-friendly sustainable case, the textured filament has expanding technology for a deep clean. Charcoal coated for adsorption and antimicrobial for healthy white teeth, the unique black coloring of the floss exposes plaque lurking in the mouth.

How do I load the BURST Floss into the BURST case?

In order to load the BURST Floss into the case you must first remove the cover over where the bobbin sits. This is done by opening the floss case and then turning the case over and pinching and pressing the 2 tabs at the same time. Turn the case back over and you will be able to lift the cover. Once this is out you can place the bobbin on the spool and thread the floss to the left and replace the cover. This sits on and clips into place.

Feel like chatting to an oral care enthusiast?

Reach out to: hello@burstoralcare.com

BURST® is a registered trademark of BURST,USA Inc.

BURST,USA Inc.
440 N BARRANCA AVE #4658 COVINA, CA 91723

www.burstoralcare.com

