Infant-to-Toddler Electric Toothbrush

EN Cleaning

Rinse brush head with water before and after each use. Inspect carefully before each use. Ensure the brush head is connected properly. Throw away at the first sign of damage or weakness. Not suitable for boiling, steam sterilization, or dishwasher use.

Instructions

Always use under adult supervision. Press and release the power button to turn on, adjust vibration strength, and to turn off. Be gentle when brushing around new teeth. Rinse the brush head with water after each use. Replace brush head every 3 months or earlier if bristles appear worn/faded or loose.

Powered by 1 AAA Alkaline, non-rechargeable battery (included). To access the battery compartment, turn the notch at the bottom of the handle in the direction of the unlock symbol using a flat edge tool. Once unlocked, remove the bottom cover. To insert a new battery, observe proper polarity as indicated. Use only AAA alkaline batteries. Dispose of old battery properly. Follow local or state regulations for battery disposal or recycling.

WARNING!

- Always use this product with adult supervision.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance. Keep out of reach of children when not in use.
- Not suitable for children under 1 year of age.
- For tooth brushing only. DO NOT BITE OR CHEW. Do not use as a toy.
- Do not submerge in water.
- Do not use if any parts are loose or damaged.
- Contains a battery. Keep battery compartment securely closed at all times.
- Exhausted batteries are to be removed from the appliance and safely disposed of.
- If the appliance is to be stored unused for a long period, the batteries should be removed.

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