

Personal Microderm Classic

Personal Microderm Classic for fresh looking, radiant and clear skin.

This revolutionary at-home exfoliation device provides the same brilliant results as professional, in-office microdermabrasion treatments. Combines patented spinning disc technology with perfectly calibrated vacuum suction to leave skin looking fresh and radiant.

Weekly treatments reduce the appearance of fine lines and wrinkles, blemishes and enlarged pores. Patented spinning discs exfoliate away dead, dull cells from the surface of the skin which stimulates new cell growth and makes the skin glowing and youthful. Your skin care products can also penetrate 20 times deeper maximizing the benefits of skin care.



Benefits

- ✓ Diminish the appearance of fine-lines & wrinkles
- ✓ Reduce the appearance of pore size
- ✓ Diminish the appearance of blemishes
- ✓ Even out skin texture
- ✓ Reduce the appearance of dark spots
- ✓ Boosts collagen and elastin production
- ✓ Stimulates new cell growth



Technical Information



How to Use



1: Cleanse and dry skin thoroughly.



2: Place disc of choice & screw on cap.



3: Plug in & turn on device.



4: Using the white disc, practice on arm or leg before using on face.



5: Move the device quickly across skin in an upward motion and keep the device constantly moving. Don't hover over any spots. Do only one pass per area per treatment.



6: Pull skin tight to help device smoothly glide over skin.



7: Follow with toner and moisturizer. Repeat every 6-7 days.



8: Wash the cap and filter after each use with water and soap.

Use only as directed. Not following instructions (including multiple passes, too high setting, pressure, pausing in one spot, or other unintended uses) may lead to rash, redness, itching or even injury. If this occurs, discontinue use, and consult a medical professional. Product is intended for external use only. Avoid sensitive areas including any area with broken, irritated, or itching skin. The product and its charger are not waterproof. Do not use near water, while skin is wet, or while skin is moist from product. Keep product out of reach of children and do not use on children below age 16. PMD is not responsible for malfunction or injury if the product is not used as intended. Not intended for clinical or commercial use. The effects of the Personal Microderm are not always immediately visible, and results may vary. Do not use any active or peeling skincare products for 48 hours before or after using the Personal Microderm — your skin may have increased sensitivity.