

Wave Pro

The PMD Wave Pro is an advanced skincare spatula that extracts, infuses, and lifts using SonicGlow™ Technology and EMS (electronical muscle stimulation) to leave the skin glowing and radiant.

PMD Wave Pro has 3 modes: Exfoliating & Extracting, Infusing and Lifting.

Extract: Extract mode uses negative ions that help lift debris, dead skin, and sebum out of the pores. Precise Xtractor™: Paired with Extract mode, Precise Xtractor™ helps to remove blackheads and whiteheads.

Infuse: Infuse mode uses positive ions that split skincare into smaller molecules. This allows skincare to be more easily absorbed into the skin.

EMS (electronical muscle stimulation): EMS uses an electrical current to stimulate the muscles, helping to give a more toned and youthful appearance.



Benefits

- ✓ Minimizes the appearance of enlarged pores and blackheads
- ✓ Reduces oil congestion
- ✓ Makes the skin softer and smoother
- ✓ Enhances the absorption of skincare products
- ✓ Lifts and tones the skin



Technical Information



How to Use



1: Dampen skin with water or toner. Turn on mode 1 (Extract). With the front side turned down, gently glide in an upward and outward motion. Place Precise Xtractor™ over blackhead or whitehead and gently push down to extract. It is safe to use Extract mode once a week.



2: Apply skincare to skin and power on mode 2 (Infuse). Using back of the device, gently infuse skincare in an upward motion. It is safe to use mode 2 daily.



3: Mist face with toner. Turn on mode 3 (Lift). Ensure your fingers are on both Touch Points on the sides of the Wave PRO and lift in an upward and outward direction by crow's feets, eyebrows, and lip lines. It is safe to use mode 3 daily.