RENPHO

USER MANUAL

SHIATSU FOOT MASSAGER

Model: RF-FM059R

CONTENTS

• ENGLISH 01~10

• DEUTSCH 11~20

• FRANÇAIS 21~30

• ITALIANO 31~40

• ESPAÑOL 41~50

Thank you for your purchase of RENPHO Foot Massager. The massager offers a powerful massage. While enjoying massager, you should use it with care and safety. To ensure safety and avoid any property damage, please read this manual carefully and use as indicated. Non-compliance may cause serious injuries or severe property damage.

IMPORTANT INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USING.

1. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
2. Children should be supervised to ensure that they do not play with the appliance.
3. **WARNING** – To reduce the risk of burns, fire, electric shock, or injury to persons.
 - a) This appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
 - b) Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
 - c) Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
 - d) Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
 - e) Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
 - f) Do not carry this appliance by supply cord or use cord as a handle.

- g) Keep the cord away from heated surfaces.
 - h) Never operate the appliance with the air openings blocked.
Keep the air openings free of lint, hair, and the like.
 - i) Never drop or insert any object into any opening.
 - j) Do not use outdoors.
 - k) Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
 - l) To disconnect, turn all controls to the off position, then remove plug from outlet.
4. The appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
 5. The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.
 6. **WARNING:** Use heated surfaces carefully. May cause serious burns. Do not use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
 7. **WARNING:** Do not stand on or in appliance. Use only while seated.
 8. The device is intended for household use only and not intended for medical use. It should not be used in treatment or alleviation of disease or abnormal physical state.

About the Battery

1. Remove and immediately recycle or dispose of used batteries according to local regulations and keep away from children.
2. Do NOT dispose of batteries in household trash or incinerate.
3. Even used batteries may cause severe injury or death.

4. Call a local poison control center for treatment information.
5. Clean the battery contacts and also those of the equipment prior to battery installation.
6. Non-rechargeable batteries are not to be recharged.
7. Do not force discharge, recharge, disassemble, heat above (manufacturer's specified temperature rating) or incinerate. Doing so may result in injury due to venting, leakage or explosion resulting in chemical burns.
8. Ensure the batteries are installed correctly according to polarity (+ and -).
9. Do not mix old and new batteries, different brands or types of batteries, such as alkaline, carbon-zinc, or rechargeable batteries.
10. Remove and immediately recycle or dispose of batteries from equipment not used for an extended period of time according to local regulations.
11. Always completely secure the battery compartment. If the battery compartment does not close securely, stop using the product, remove the batteries, and keep them away from children.
12. Do not expose the appliance or battery to excessive temperatures.
13. Always select the correct size and type of battery most suitable for the intended use. Information provided with the equipment to assist correct battery selection should be retained for reference.
14. In case of ingestion of a cell or battery, seek medical assistance promptly. Swallowing lithium coin cells or batteries can cause chemical burns, perforation of soft tissue, and in severe cases can cause death. They must be removed immediately if swallowed.
15. Do not allow children to replace batteries without adult supervision.
16. Do not force discharge batteries. Do not deform batteries. Do not dispose of batteries in fire.

17. Exhausted batteries should be immediately removed from equipment and properly disposed of. When discharged batteries are kept in the equipment for a long time, electrolyte leakage can occur causing damage to the equipment and/or personal injury.
18. A lithium battery with a damaged container should not be exposed to water.
19. Store unused batteries in their original packaging away from metal objects. If already unpacked, do not mix or jumble batteries.
20. Do not use modified or damaged batteries.
21. Be aware of the risk of terminals of the battery-operated appliance or battery being short- circuited by metal objects.

SAVE THESE INSTRUCTIONS

About the RENPHO Foot Massager

1. Specifications

Model: RF-FM059R

Input: 12V= 2.5A

Rated power: 30W

Heating area: Bottom

Time: 15min/30min

Material: ABS + Lycra fabric

Heating temperature: Blow 50 /122

Color: Black / Silver / White

2. What's in the box

Shiatsu Foot Massager x 1

Remote Control x 1

Adaptor x 1

3. Product Introduction

4. Function

Time: Select 15 or 30 minutes automatic shut-off.

Knead Intensity: Select the kneading intensity among low (L), medium (M), high (H) or OFF.

Air Intensity: Select the air pressure intensity among low (L), medium (M), high (H) or OFF.

Heat: Press to turn on/off the heating function.

Using the RENPHO Foot Massager

Enjoying a foot massage for 15 or 30 minutes once or twice daily can relieve muscle pain, improve blood circulation, and effectively improve the symptoms of insomnia. And we recommend you limit massaging to 15-30 minutes. Longer massage times may cause mild pain.

1. Connect the plug to the power source. The Power indicator should flash and the massager is in standby.

2. Please remove the protective film from the touch panel before use.

Press the Power button. The knead and air intensity will be at default low level (L). The default auto shut-off will be 15 minutes with the indicator light red.

3. Press the Time button to change the auto shut-off to 30 minutes and the indicator light will turn blue.

4. Press the Knead Intensity and Air Intensity buttons to adjust the knead and air pressure intensity (L, M, H, OFF).

5. Press the Heat button to turn on/off the heating function.

Note: The heating function cannot work alone. The air pressure and kneading functions can work independently. The heating function can only work when the air pressure and/or kneading massage are working.

TIPS: Use your thumb to press buttons to get better touch sensitivity.

6. The foot massager fits up to US Men's size 12 (EU Men's size 46).

7. To remove and wash the foot cover:

1

2

3

4

Product Care and Maintenance

1. The use of non-rated voltage not (DC12V/2.5A) is prohibited.
2. It is prohibited for non-professional personnel to repair the massager or replace parts.
3. If the power adapter is damaged, it must be replaced by the maintenance department or similarly qualified persons approved by the manufacturer so as to avoid potential hazards.
4. Do not use or store it in wet environments or high temperatures.
5. Unplug the power adapter from the outlet after using.
6. Ensure that the operation of this product does not exceed more than 30 minutes at any given time.
7. If you feel unwell or in any unusual discomfort, please stop using the product immediately.
8. If high temperature or foul odor is experienced while using the product, stop using immediately and unplug the device.
9. Any Individual who may be pregnant or recently gave birth, has a heart condition or a pacemaker, suffers from diabetes or other sensory impairment, phlebitis and/or thrombosis, joint dysfunction, inflammation or injuries, weak bones, osteoporosis, a high fever, an abnormal or curved spine, is at an increased risk of developing blood clots, has pins/screws/artificial joints or other medical devices implanted in his/her body, is being treated for a medical condition or physical injury or is under a doctor's care should consult with a medical professional before using a massaging device designed for home use.
10. Do not use sharp objects to handle the product in order to avoid damage.
11. Please keep the product clean. You can use a soft cotton cloth and neutral detergent to clean it.
12. Do not stand when you are using the massager. Do not place excessively heavy objects on the device in order to avoid damage.

Troubleshooting

Common problems	Cause	Solutions
The product does not work.	1. Disconnected power or loose plug causing poor contact. 2. The power switch is off or the timer has been reset. 3. The motor auto shuts-off due to the device being on too long activating overheating protection.	1. Re-plug to the power source. 2. Turn the power switch on again. 3. Use the massager after the motor cools off naturally.
The massage intensity is too high or too low.	1. Your foot is positioned at the wrong angle. 2. Didn't choose intensity level which right for you.	1. Adjust the contact angle to get the moderate intensity. 2. The massage has 3 intensity levels, choose the mode you prefer. We recommend low mode to begin with.
No heating or insufficient heating.	The heat function is manufactured according to International safety standards, which limits how fast the product can heat up. Therefore, the heating-up may seem a bit slow, 7-10 minutes before you can feel the obvious heat.	If you press the foot cover to see the light, the heating is turned normally. The rate of heating is related to environment. If you are in a cold area, the heating time may take longer.
The touch panel is not working properly.	1. The plastic film on the touch panel is not removed. 2. Fingers are wet. 3. Insufficient touch area.	1. Remove the plastic film before use. 2. Dry your fingers before touching the button. 3. Use your thumb to press buttons to get better touch sensitivity.

Warranty Policy

Your RENPHO product purchase is covered by a one-year limited manufacturer warranty from the date of delivery.

For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty. If you choose not to register your product, it will not diminish the product warranty.

Customer Service

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

 **Email:** support-eu@renpho.com (UK&EU)
support@renpho.com (US&CA)
support-au@renpho.com (AU)

 **TEL:** **+44(0)7434 666088 (UK ONLY)**
Monday-Friday 9:00AM-5:00PM (GMT)
+1(844) 417 0149 (US&CA)
Monday-Friday 9:00AM-4:30PM
1800 MY RENPHO (1800 69 7367) (AU ONLY)
Monday-Friday 9:00AM-6:00PM (AEST)

*For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

1. Meaning of the crossed-out wheeled bin symbol

The crossed-out wheeled bin symbol means that you are legally obliged to dispose of these appliances separately from unsorted municipal waste, i.e. old appliances must be disposed of separately from household waste. Batteries or rechargeable batteries that are not permanently installed must be removed beforehand and disposed of separately.

2. Options for returning old appliances

Owners of WEEE can return or collect it within the framework of the possibilities set up and available by public waste management authorities to ensure proper disposal.

1. The symbol of the crossed-out garbage can on batteries or accumulators indicates that they must not be disposed of in household waste at the end of their life.
2. You are legally obliged to return old batteries and accumulators after use. You can do this free of charge at a retail store or other collection point in your area. Addresses of suitable collection points can be obtained from your city or local government.
3. Batteries may contain substances that are harmful to the environment and human health. The separate collection and recycling of old batteries and accumulators is intended to avoid negative effects on the environment and human health.
4. Please avoid the generation of waste from old batteries as far as possible, e.g. Please avoid littering public spaces by not carelessly leaving batteries or electrical and electronic equipment containing batteries.

Hereby, Shenzhen Ruiyi Business Technology Co., Ltd. declares that this equipment complies with Directive 2014/30/EU, Directive 2014/35/EU, Directive 2011/65/EU, The Electrical Equipment (Safety) Regulations 2016, The Electromagnetic Compatibility Regulations 2016, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, and other relevant provisions. The full text of the Declaration of Conformity can be obtained at the following internet address:
<https://renpho.eu/pages/compliance-data>

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FCC Regulatory Compliance

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED Regulatory Compliance

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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