

ThermoWorks WAND™ Mini

Quick Reference

Measuring Forehead Temperatures Quick, Easy, and Reliable!

- 1 Power on device**
- 2 Point at center of forehead 0.5-1.5 inches away**
- 3 Press START for 1 second—two short beeps tells you it's done**



WAND Mini is pre-set to display temperature in Fahrenheit. Three short beeps followed by one long beep while taking a forehead reading means that the temperature exceeded 99.5°F (37.5°C). If the temperature reads above 100.4°F (38°C), the red LED will light up. Please wait 2 seconds between measurements.



ThermoWorks

Technical Support

For warranty, service, and technical assistance, please contact ThermoWorks' Technical Support at (385) 330-0591 or email techsupport@thermoworks.com.

Functions

Power Off Device – Device will auto-off after 60 seconds to save battery life. Or, press and hold Power button for about 10 seconds. 'Off' will flash before powering off.

View Readings – With device turned on, press Start button momentarily to see the last reading.

Change from °F to °C – With instrument turned off, press and hold the START button. After about 8 seconds, the display will show "°F" flashing. Release START button and press again and release. The display will show "°C" flashing. After about 5 seconds, the display will auto reset and there will be 2 short beeps. Complete the process again to switch back to °F. Please note: stored readings will be deleted when switching between units.

Cautions

- Forehead should be clean, dry, and unobstructed before taking readings.
- WAND Mini is not waterproof. Avoid wet areas and splashes.
- Do not expose device to temperatures above 100°F (38°C) or below 50°F (10°C).
- Replace batteries when battery icon flashes.
- Store device at room temperature where readings will be taken for best results.
- Avoid drops, bumps, and rough handling that could damage device.
- Sensor lens area is the most delicate part; keep clean and free from debris.
- Consult user manual before cleaning lens.