medela

Instructions for use

In-bra Breast Milk Collector

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In-bra Milchauffangschalen

Mode d'emploi

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In-bra

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In-bra

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Bh-moedermelkcollector

Bruksanvisning

In-bra Bröstmjölksuppsamlare

Bruksanvisning

In-bra Brystmelkoppssamler

Brugsanvisning

In-bra Brystmælksopsamler

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- Disinfection**
 - Clean the parts according to the cleaning instructions.
 - Cover the parts with water and boil for 5 minutes in a pot.
 - Alternatively, the In-Bra Breast Milk Collector can be cleaned in the dishwasher on the top rack only. It is safe to disinfect in the Medela Quick Clean™ Microwave Bags.
 - Place the separated parts on a clean dry cloth to air-dry.

Material Shell (front-side): Polypropylene
Material Membrane (back-side): Silicone

Disposal
According to local regulations.

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ENGLISH

The Medela In-bra Breast Milk Collector collects every precious drop of milk and ensures a leak-free experience thanks to its plug. The anatomic design creates a gentle feel on breast and makes them comfortable and discreet to wear.

Indication
- To be worn during breastfeeding or pumping, when the milk flow starts in the other breast - To be worn when too much milk leaks between breastfeeding or pumping session
Assembly and Use
- Position the arrow on top of the silicone membrane to the notch on top of the shell. Connect the silicone membrane onto the shell. Be sure it snaps on all the way around. - Position the assembled Milk Collector inside your bra with the opening centered over your nipple. Ensure the pouring spout is facing upwards. We recommend wearing a nursing bra, as the stretch fabric ensures less pressure on the breast. - When moving around whilst wearing the Milk Collector, close the pouring spout with the plug on the top to avoid spillage. - To remove the Milk Collector, slightly lean forwards, grasp the groove and remove upwards. Place the Milk Collector flat base on the table surface to avoid spillage. - To store the collected drip milk, remove the plug from the pouring spout and fix it into the hole on the silicone band. Pour the milk to a storage container such as Medela Easy Pour Breast Milk Storage Bag or bottle. For information on storing your breast milk, visit https://www.medela.com/storing-breast-milk. - Condensation may appear in the collector. When at rest, to allow ventilation, open the spout and fix the plug into the hole on the silicone band.
Caution: Not to be worn when sleeping or in a reclined flat position. Replace the In-bra Breast Milk Collector with a nursing pad when sleeping. When used for an extended period, remove the milk collectors every 2-3 hours for feeding/pumping to allow air to circulate and avoid breast compression. Milk collected between breastfeeding/pumping should be discarded.
Cleaning Instructions
Before first use and subsequently once daily - Clean the parts (silicone membrane and shell) as described under "After each use". - Disinfect the parts as described under "Disinfection".
After each use - Always wash all parts as soon as possible after each use. - Rinse milky parts with cold running drinking quality water. Wash the parts in warm water with mild dish soap. Rinse parts with cold water. - Place the separated parts on a clean dry cloth to air-dry.

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