
RENPHO

User Manual

U-NECK MINI

Model: R-L001

Contents

· ENGLISH	01 ~ 11
-----------	---------

· DEUTSCH	12 ~ 21
-----------	---------

· FRANÇAIS	22 ~ 31
------------	---------

· ITALIANO	32 ~ 41
------------	---------

· ESPAÑOL	42 ~ 51
-----------	---------

IMPORTANT SAFETY INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USING.

1. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
2. Children should be supervised to ensure that they do not play with the appliance.
3. The device is intended for household use only and not intended for medical use. It should not be used in treatment or alleviation of disease or abnormal physical state.
4. Please consult a doctor before use if you fit into any of the following user groups:
 - Persons who are pregnant or postpartum
 - Persons with skin infections, allergies, redness and/or skin burns
 - Persons with osteoporosis
 - Persons who had undergone surgery
 - Persons with any swelling or wound around the abdomen
 - Persons with a pacemaker or other implantable medical devices
 - Patients with heart disease, malignant tumors, or other serious health conditions
 - Persons who are receiving medical treatment
 - Persons with heart problems, abnormal blood pressure, diabetes or neurological disorders
5. If you feel unwell or in any unusual discomfort, please stop using the appliance immediately.
6. NOT suitable for children.
7. DO NOT use it while driving.
8. The elderly need to use under the doctor's guidance.
9. DO NOT use this appliance after drinking or taking sleeping pills.
10. NEVER try to inspect or open the appliance in any way.
11. DO NOT use this appliance for too long, as this may cause minimal damage to the skin.
12. DO NOT put heavy objects on the appliance.

13. NEVER use pins or other metallic fasteners with this appliance.
14. DO NOT use this appliance if you apply any oils or lotions to your skin.
15. NEVER clean this appliance unless it is unplugged. Use a cloth dampened with a very small amount of water to clean. You can use non-abrasive cleaners. DO NOT plug this appliance back in until it is totally dry.

CAUTION - To reduce the risk of burns:

1. The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.
2. Use heated surfaces carefully. May cause serious burns. Do not use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.

WARNING - To reduce the risk of burns, fire, electric shock, or injury to persons:

1. This appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
2. Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
4. Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
5. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
6. Do not carry this appliance by supply cord or use cord as a handle.
7. Keep the cord away from heated surfaces.
8. Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, and the like.

9. Never drop or insert any object into any opening.
10. Do not use outdoors.
11. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
12. To disconnect, turn all controls to the off position, then remove plug from outlet.
13. The appliance is only to be used with the power supply unit provided with the appliance.
14. The appliance must only be supplied at safety extra low voltage corresponding to the marking on the appliance.

SAVE THESE INSTRUCTIONS

WHAT'S IN THE BOX

U-Neck Mini x1

USB-C Charging Cable x1

User Manual x1

SPECIFICATIONS

Product Name: Neck Massager


Model: R-L001

Input: 5V = 2A

Rated Power: 10W

Battery Capacity: 2000mAh 14.8Wh

ABOUT YOUR U-NECK MINI

*During use, gently pull the positioning handles downward for an optimal fit, ensuring the massager nodes are in close contact with the desired massager areas.


INDICATOR

Status	Indicator
Device turned on / Wireless Connectivity disconnected	Turns solid green
Wireless Connectivity connected	Turns solid blue
Low battery / Charging	Flashes red slowly
Fully charged	Turns solid red
Charging at a voltage higher than 5V	Keeps flashing red quickly for 1 minute and then turns off, indicating the wrong charger is used.

CHARGING YOUR U-NECK MINI

Fully charge the device before the first use.

Charging Time: approx. 3 hours

Note: Do not use the device while charging.

USING YOUR U-NECK MINI

Tip: Tie up long hair, remove any neck accessories and make sure the massage node covers are securely attached before wearing the device. Ensure comfortable wearing before turning on the device.

Button/Feature	Description
 Power Button	Press and hold for 2 seconds to turn on/off the device. After power on, low heat, weak intensity, built-in music and Wireless Connectivity are enabled by default.
 Intensity Button	Press to switch between weak intensity and strong intensity.
 Heat Button	Press to adjust the heat level between low heat and high heat or to turn off the heating function.
 Music Button	a. Press to switch to the next track. b. Double press to turn on/off the voice broadcast. c. Press and hold to turn on/off the music.
10-Minute Auto-off Timer	The device will automatically shut off after 10 minutes of continuous use. Restart the device if needed.

CONNECTING TO WIRELESS CONNECTIVITY

1. Turn on the device and enable Wireless Connectivity on your mobile device.
2. Search for “**R-L001**” via your Wireless Connectivity settings and connect to play your music playlist.

Note: Once paired successfully, the device will be automatically connected to Wireless Connectivity next time, followed by a voice prompt “Wireless Connectivity is connected”.

CLEANING & MAINTENANCE

1. Turn off the device and allow it to cool down before cleaning or storing.
2. Wipe the device surface with a soft, dry cotton cloth. If there are stains or stubborn dirt, use a cotton ball or cloth soaked in alcohol to clean. Do not submerge the device in water or use bleach or corrosive cleaners.
3. To clean the massage node covers, gently detach them from the massage nodes, hand wash with mild detergent and rinse thoroughly with clean water.
4. Allow the device and massage node covers to air dry completely before use or storage.
5. Store the device in a cool, dry place. Do not expose it to direct sunlight, moisture, heat, and corrosive gases.
6. If you plan to store your device for longer than six months, charge it to 50% every six months.
7. Do not unzip, disassemble or repair the device by yourself as this may cause damage and void the manufacturer's warranty.

TROUBLESHOOTING GUIDE

Problem	Possible Cause	Solution
Power-on failure.	Short press the Power Button.	Press and hold the Power Button for 2 seconds to turn on the device.
	Out of power.	Charge the device.
Charging failure.	Use other types of charging cables.	Use the provided charging cable only.
	Use a charger with voltage higher than 5V.	Use the DC 5V 2A charger.
The device shuts off suddenly.	Out of power.	Charge the device.
	10 minutes auto-off timer is activated.	Restart the device if needed.
No voice broadcast.	The voice broadcast is turned off.	Double press the Music Button to turn on the voice broadcast.
No music.	The music is turned off.	Press and hold the Music Button to turn on the music.
	The Wireless Connectivity music volume is at level 0.	Adjust the Wireless Connectivity music volume via your mobile device.
If you encounter any other problems, please contact our customer service.		

WARRANTY POLICY

Your RENPHO product purchase is covered by a one-year limited manufacturer warranty from the date of delivery.

For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty.

If you choose not to register your product, it will not diminish the product warranty.

CUSTOMER SERVICE

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

 TEL: +44(0)7434 666088 (UK ONLY)
Monday-Friday 9:00AM-5:00PM (GMT)
+1(844) 417 0149 (US&CA)
Monday-Friday 9:00AM-4:30PM
1800 MY RENPHO (1800 69 7367) (AU ONLY)
Monday-Friday 9:00AM-6:00PM (AEST)
+49 40 23969285 (DE ONLY)
Monday-Friday 9:00AM-6:00PM (CET)

 Email: support-eu@renpho.com (UK&EU)
support@renpho.com (US&CA)
support-au@renpho.com (AU)

*For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

FCC Regulatory Compliance

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED Regulatory Compliance

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1: This device may not cause interference.

2: This device must accept any interference, including interference that may cause undesired operation of the device.

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1. Meaning of the crossed-out wheeled bin symbol

The crossed-out wheeled bin symbol means that you are legally obliged to dispose of these appliances separately from unsorted municipal waste, i.e. old appliances must be disposed of separately from household waste. Batteries or rechargeable batteries that are not permanently installed must be removed beforehand and disposed of separately.

2. Options for returning old appliances

Owners of WEEE can return or collect it within the framework of the possibilities set up and available by public waste management authorities to ensure proper disposal.

3. Personal data

All end-users of WEEE are reminded that you are responsible for deleting personal data from the WEEE to be disposed of.

1. The symbol of the crossed-out garbage can on batteries or accumulators indicates that they must not be disposed of in household waste at the end of their life.
2. You are legally obliged to return old batteries and accumulators after use. You can do this free of charge at a retail store or other collection point in your area. Addresses of suitable collection points can be obtained from your city or local government.
3. Batteries may contain substances that are harmful to the environment and human health. The separate collection and recycling of old batteries and accumulators is intended to avoid negative effects on the environment and human health.
4. Please avoid the generation of waste from old batteries as far as possible, e.g. Please avoid littering public spaces by not carelessly leaving batteries or electrical and electronic equipment containing batteries.

Hereby, Shenzhen Ruiyi Business Technology Co., Ltd. declares that this equipment complies with Directive 2014/53/EU, Directive 2011/65/EU, Radio Equipment Regulations 2017, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, and other relevant provisions. The full text of the Declaration of Conformity can be obtained at the following internet address: <https://renpho.eu/pages/compliance-data>

Maximum RF power: 0.65dBm
Operation frequency: 2402MHz to 2480MHz

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