tria®

Tria 4X

ADVANCED LASER TECHNOLOGY


User guide & manual.

Skin tone chart

Check your skin tone to find out if you can safely use the Tria 4X. Match the skin in the area you wish to treat to the range of skin tones shown on this chart. If your skin tone matches any of the skin tones labelled Not Safe for treatment, then DO NOT USE the Tria 4X.

Intended skin tones:

Not safe:

FOR USE ON: Legs; Arms; Bikini Line; Underarms; Feet; Hands; Chest; Back; Stomach; Women's Facial Hair (cheek line downward)

Read these instructions before use

Before you use your new Tria 4X, read this booklet for important information:

- Confirm that the Tria 4X is safe for you.
- Learn how to use the Tria 4X.
- View our treatment results timeline.
- Read our tips for best results.
- Review contraindications and warnings, which are marked with the symbol below:

Be sure to read and follow all the warnings and safety information. We recommend that you save these instructions for future reference. If you have questions or concerns, please contact Customer Service.

Contents

Introduction	4
What's in the box?.....	5
Indications for use.....	6
Contraindications	7
Warnings	8
Charging.....	12
How to Use	14
What to Expect	15
Device Care	16
Troubleshooting.....	17
Warranty.....	20
Technical Specifications	21
Français	25
Deutsch	45
Español.....	65
Italiano	85
繁體中文	105
简体中文	125
日本語	145
한국어.....	165

Introduction

Pioneers in laser-based skincare for over 20 years, Tria was founded by scientists who set the 'gold standard' for in-clinic laser hair removal. Their mission was to empower you to achieve those same clinical-grade results with a handheld laser from the comfort and privacy of your own home.

Our priority has always been to give you excellent results, utilising the best technology. Our multi-award-winning laser devices utilise the same diode laser technology that is preferred by professionals around the world, so you can be confident, knowing that your Tria device is as close to in-clinic technology as possible in a handheld device.

Get clinically proven, permanent results* at home with your Tria 4X. Subjects reported up to 100% hair removal in high fluence treatment areas following a full treatment regimen.** Your Tria 4X expertly targets the melanin in your hair, disabling hair follicles and preventing future hair growth.

*Permanent results are defined as 'a permanent reduction in hair regrowth'. i.e., a long-term stable reduction in hair counts following full course of laser treatments.

**Based on subject evaluations reported one year post-treatment, 31% of subjects reported 100% hair removal at the highest power setting during an independent clinical trial involving 21 subjects.

What's in the box?

TRIA 4X
Power Supply
IFU
Quick Start Guide (QSG)

Indications for use

The Tria 4X is used to remove unwanted hair. It is an over-the-counter system intended for adjunctive use with shaving for hair removal sustained with periodic treatments. It is also intended for permanent reduction in hair regrowth defined as long-term, stable reduction in the number of hairs regrowing when measured out to 6, 9 and 12 months after the completion of the treatment regime. It is intended for use only on light to medium skin tones (see the Skin Tone Chart on the inside front cover) and naturally brown or black hair. It will not work on white, grey, blonde or red hair.

Contraindications

DO NOT use the Tria 4X if your skin tone matches one listed as Not Safe on the Skin Tone Chart on the inside front cover of this booklet. If you use the Tria 4X on skin deemed Not Safe for treatment, the skin may absorb too much laser light and be injured (burned, discoloured or scarred).

Please Note: The skin sensor is a safety feature that unlocks the Tria 4X. If it is not safe for you to use the Tria 4X, you won't be able to unlock it with the skin sensor.

DO NOT use the Tria 4X if you are under the age of 18 or use it on anyone else under the age of 18. The use of the Tria 4X on or by children has not been studied and could result in serious injury.

DO NOT use the Tria 4X if you have a lighter hair colour, such as blonde, white, red or grey. In order for the laser to work, you must have naturally dark hair where you are treating. Lighter hair pigments will not absorb enough of the laser's energy to disable the hair follicle.

Warnings

DO NOT use the Tria 4X on or around the eyes, eyebrows or eyelashes. Doing so can cause permanent eye injury.

DO NOT use the Tria 4X if the Skin Tone Chart shows that your skin tone is not suitable. If you have African, Indian, Native American or Pacific Island heritage, you may have a skin tone that is not suitable — be sure to check the Skin Tone Chart. Darker skin tones can absorb too much laser light, which can injure the skin.

DO NOT use the Tria 4X anywhere on the ears, nipples, genitals or around the anus. These areas may have more sensitive skin, a darker skin tone, and/or a greater hair density, and using the Tria 4X there may cause serious skin injury.

DO NOT use the Tria 4X on any area where you may someday want your hair to grow back. The results may be permanent and irreversible.

DO NOT use the Tria 4X on parts of your skin that are darker than your usual skin tone: tattoos, dark brown or black spots (such as large freckles, birthmarks or moles) or skin that is red or swollen from a previous treatment. Dark skin may absorb too much laser light, which can injure the skin.

Hair removal by laser can cause increased hair growth in some individuals. Based upon currently available data, the highest risk groups for this response are females of Mediterranean, Middle Eastern and South Asian heritage treated on the face and neck.

DO NOT use the Tria 4X on damaged skin, including open wounds, cold sores, sunburns or areas where you have had a skin peel or other skin resurfacing procedure within the last 6-8 weeks (check with your doctor). These conditions can make the skin more susceptible to injury.

DO NOT use the Tria 4X if you are pregnant or breast-feeding. The Tria 4X has not been tested on pregnant women. Hormonal changes could increase skin sensitivity and the risk of skin injury.

DO NOT use the Tria 4X if it appears damaged, and DO NOT attempt to open or repair the Tria 4X. These actions increase the risk of skin or eye injury. Please contact Customer Service for assistance.

Always follow the basic safety precautions for electrical devices:

DO NOT get the Tria 4X or the charger damp or wet. DO NOT use the Tria 4X while bathing or put it next to a bath or sink. DO NOT operate a charger with a damaged plug or cord. ALWAYS keep the cord away from heated surfaces. Failure to follow these precautions can cause electrocution or electric shock.

DO NOT use the Tria 4X for any purpose other than hair removal.

The Tria 4X is for personal use. The Tria 4X is not for use in a medical office, clinic, salon or other commercial setting.

Invisible Laser Radiation. Do Not View Directly with Optical Instruments.
Class 1 Laser Product.

Viewing the laser output with certain optical instruments (for example eye loupes, magnifiers, and microscopes) within a distance of 100 mm may pose an eye hazard.

No modification of this equipment is allowed.

This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Benefits and risks

A medically supervised clinical study has shown that the Tria 4X is safe and works when used as directed. Seventy-seven (77) subjects self-administered a 6-week course of treatment (3 treatments). Subjects were allowed to choose their desired fluence (7, 12, 20J/cm²) and all side effects were observed. The number of hairs at each treatment site was counted before each treatment and at 1, 2, 3, 6 and 12 months after the last treatment. No additional treatments were performed after the third treatment.

Benefits

When used as directed, the Tria 4X gives you permanent results. The results vary from person-to-person and depend on which body area is treated and how closely treatment instructions are followed.

The clinical study demonstrated the following benefits:

- 61% average reduction in hair count 3 weeks after the first treatment
- 70% average reduction in hair count 3 weeks after the second treatment
- 60% average reduction in hair count 1 month after the third treatment
- 33% average reduction in hair count 12 months following the third treatment

In the study, 87% of subjects reported finer hair at 1 month after the last treatment than before any treatment. And 70% of subjects reported lighter hair at 1 month after the last treatment than before any treatment. These results were obtained with only 3 treatments. For best results, treat every 2 weeks for 3 months or until you are satisfied with the amount of hair reduction.

Risks

It is important to use the Tria 4X properly. Make sure that you understand and follow all warnings (see pages 5-7). The skin sensor will prevent use of the device on contraindicated skin types (see page 5).

If you use the Tria 4X as directed, the risks are low:

- In the clinical study, almost all users reported mild to moderate pain or discomfort during their first treatment, with a reduction in sensitivity reported in their later treatments. This is normal. You may feel warmth, burning, tingling or itching; some users described feeling the same level of discomfort as when a rubber band snaps against the skin. If the pain is too intense or persists after a treatment, stop using the Tria 4X and check with your doctor before using it again.
- Your skin may become red right after treatment or within 24 hours. This redness generally clears up within a few hours. See your doctor if the redness does not go away within 2–3 days.
- In rare cases, some users may see an increase in the number of hairs in the treated area after the initial treatment. If this increase persists over several treatments, stop using the Tria 4X and see your doctor.

The following symptoms affected no users in the clinical study (0 out of 88) but have been reported with other hair removal lasers:

- Skin swelling, especially around the hair follicle. If the swelling doesn't go away within 2–3 days, stop using the Tria 4X and see your doctor.
- Temporary or permanent skin discolouration. If your skin colour changes, stop using the Tria 4X and see your doctor.
- Blisters. This is very unlikely for people with a light or medium skin tone. If you do develop a blister, stop using the Tria 4X and don't let the blister become infected. Consult your doctor if needed.

Charging your device

1. The device may come with some charge, but we always recommend fully charging before first use.
 2. Insert the cord into the port on the device. Plug the power supply into an outlet.
 3. When your device is charging, the battery indicator light will roll.
 4. When fully charged, the indicator will remain solid.
- NOTE: To prevent damage to your Tria 4X, use only the included Tria battery charger. If you lose the charger, please contact Customer Service.
 - NOTE: You may charge the laser in Quiet mode (charges without the sound of cooling fan), but it will take longer. To change charging modes, press and hold the power button while charging until you hear a tone.

Sensitivity test

1. Ensure your skin is clean and dry in the area you wish to treat.
2. Turn on your device by single pressing the power button.
3. Place the treatment window flat against your skin and slowly glide side-to-side for several seconds.
4. If Level 1 was tolerable for you, then move to another test area and increase to Level 2, testing a patch at that level with one laser pulse. To adjust the Tria 4X setting up or down, press the power button.
5. Wait 24 hours and examine the tested area. If the skin appears normal, proceed with treatment.

Retest your skin if your skin tone changes dramatically (due to tanning).

If you notice blistering, burns, changes in your skin tone or any other side effects, or if you are not sure about the outcome of the test, DO NOT USE your Tria 4X. Contact Customer Service or your doctor.

How to use

1. Ensure your skin is clean and dry in the area you wish to treat.
2. Shave the treatment area.
3. Single press the power button to turn on the device. Your device will be locked.

To unlock:

- Hold the skin sensor located on the base of your device against the area you want to treat.
 - If your skin tone is compatible, the device will beep.
 - If your skin tone is not compatible, the device will remain locked.
4. Select your treatment level from 1-5 by single pressing the power button. We recommend using the highest level you can tolerate.

5. Hold the treatment window against your skin until you hear two beeps.
- If the device does not beep, make sure that the window is flat against the treatment area.
 - If the device beeps and then buzzes, the laser pulse has been disrupted. Repeat this pulse.

6. Lift the laser and place back down in an interlocking pattern.

DO NOT repeatedly treat the same spot. Always move the treatment window to a new area after each pulse.

7. The pulse counter will count up with each successful laser pulse. Some areas may require more pulses than others.
8. Turn off the device by holding down the power button.

Repeat treatment every 2 weeks for up to 3 months or until you reach the desired results.

NOTE: The cooling fan may go on and off during your treatment. This is normal and nothing to worry about.

What to expect

During use

- The sensation varies, depending on the person, treatment area and treatment level. Some users may feel the sensation of a rubber band snapping on the skin. Others will feel a warm, prickling sensation.
- As the amount of hair in the treated areas decreases, the sensation will also decrease. For this reason, you may be able to use a higher treatment level after your initial treatments.
- Some users experience mild redness that goes away within 24 hours.

During first 3 months

- After your first treatment, there will be no visible change to the hair. Laser hair removal does not cause hair to immediately fall out.
- With additional treatments, there will be a gradual reduction in the number of hairs with each treatment.
- The hair that does regrow will also begin to lighten and become finer.

After 3 months

- Hair follicles that have been deactivated by the Tria 4X will not grow back.
- Do occasional touch-ups or additional treatments if necessary.

Device care

The Tria 4X doesn't require special care, but we do recommend that you keep it clean, especially the treatment window and skin sensor.

Cleaning

- Make sure that the Tria 4X is turned off and unplug the battery charger from the Tria 4X and the power outlet.
- To clean the Tria 4X and battery charger, moisten a soft cloth with water and mild soap, and wipe the treatment window and the exterior of the device. Be careful not to damage the treatment window. Allow to dry completely before charging or turning on your Tria 4X.
- To clean the skin sensor located at the base of the Tria 4X, wipe it gently with a dry, lint-free cloth.
- DO NOT get the Tria 4X or battery charger damp or wet. If it gets wet, DO NOT use the Tria 4X or the battery charger. Contact Customer Service. Performing a treatment with a wet Tria 4X (or component) can cause electrocution or electric shock.

Recycling/disposal

The Tria 4X contains lithium batteries that cannot be replaced. To recycle or dispose, please follow the regulations where you live or contact Customer Service for information.

Trouble-shooting

Please contact Customer Service for additional help.

My Tria 4X won't turn on.

- Your Tria 4X is not charged.
- Charge your Tria 4X.

While I'm charging my Tria 4X, the battery indicator on the display doesn't fill up.

- Your Tria 4X is not connected properly.
- Confirm that your Tria 4X is firmly plugged into the charger and the wall outlet.
- Your Tria 4X is already fully charged and ready to use.
- Unplug the battery charger from the wall and from the Tria 4X.

Even after I charge my Tria 4X, the battery indicator is still filling up.

- Your Tria 4X isn't fully charged.
- Charge your Tria 4X.

When I put my Tria 4X against my skin, nothing happens.

- Your Tria 4X is locked.
- The Tria 4X will lock when it's not in use. Unlock your Tria 4X with the skin sensor.
- Your Tria 4X is plugged into the battery charger.
- The Tria 4X will not work when it is plugged in to the battery charger. Unplug your Tria 4X.
- The treatment window isn't completely in contact with your skin.
- Confirm that the treatment window is completely flat against your skin.
- The area you're trying to treat is too small, curved or bony to treat.
- Try again, pulling the skin flat. If you can't get the whole treatment window to touch the treatment area, then the Tria 4X won't work.

I can't use my Tria 4X and an image of a lock shows on the display.

- Your Tria 4X is locked.
- The Tria 4X will lock when it's not in use. Unlock your Tria 4X with the skin sensor.

The display is flashing quickly.

- Your Tria 4X is not working properly.
- Please contact Customer Service for assistance.

Dirt or debris is stuck to the treatment window.

- Your Tria 4X needs to be cleaned.
- See cleaning instructions on page 14. If you can't get the area clean, contact Customer Service for assistance.

I tried to unlock the Tria 4X, but it's still locked.

- It may not be safe for you to use the Tria 4X.
- Match your skin tone to the Skin Tone Chart on the inside cover of this booklet, and contact Customer Service with any questions.
- You didn't use the skin sensor correctly.
- Try again, holding the skin sensor flat and firm against your skin.

When I reposition the Tria 4X, it seems to take a long time before the next pulse of laser light.

- The room is warm. The Tria 4X may pulse at a slower rate in a warm environment.
- Use the Tria 4X in a cooler environment.

My Tria 4X was working properly, but then it stopped working. The display is off and does not turn on even when I press the power button.

- The battery needs to be charged.
- Charge your Tria 4X.

When I'm using my Tria 4X, I can hear a fan working.

- This is normal.
- Continue your treatment.

Warranty

Under normal use this product is guaranteed for 2 years after its original purchase date against manufacturing defects in materials or workmanship. This warranty extends only to the original purchaser. If the product should become faulty within the warranty period, contact The Beauty Tech Group.

- The Beauty Tech Group reserves the right to repair a defective product.
- The Beauty Tech Group reserves the right to replace a defective product with the most compatible product currently available.

This warranty does not cover damage caused by the following:

- Accident, misuse, abuse or alteration
- Servicing by unauthorised persons
- Connection to incorrect power sources
- Use with unauthorised accessories

The Beauty Tech Group shall not be responsible for any incidental, special or consequential damages resulting from the use of this product. All implied warranties, including but not limited to implied warranties of fitness and merchantability are limited in duration to two years from date of original purchase. This warranty gives you specific legal rights, and you may also have other rights that exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so some of the above limitations and exclusions may not apply to you. This warranty is given in lieu of all other warranties, written or oral, whether expressed by affirmation, promise, description, drawing, model or sample. Any and all warranties other than this one, whether expressed or implied, including implied warranties of merchantability and fitness for a particular purpose are hereby disclaimed. If the product is defective within the warranty period please contact The Beauty Tech Group.

Technical specifications

Laser type:	Diode laser
Laser class:	1 According with IEC 60825-1:2014
Output wavelength:	810 nm
Output fluence:	7-20 J/cm2
Battery charger input:	100-240 V ~ , 50/60Hz, 29 VA
Battery charger output:	5.0 V DC, 2.0 A
Expected service life:	3 Years minimum

Physical

Physical operating temperature:	↕ 5°C to 30°C
Operating relative humidity:	Non-condensing 15-85%
Operating atmospheric pressure:	700 hPa to 1060 hPa
Transport and storage temperatures:	↕ -20°C to 50°C
Transport relative humidity:	Non-condensing 0-93%

Symbols

Manufacturer

“WEEE (Waste Electrical and Electronic Equipment)”. The waste products should be hand RGB legally.

Refer to instruction manual/booklet.

Protection against ingress of solids, it means the volume of the solids should less than 12mm. Protection against ingress of water, it means the device could protected against vertical falling water drop 15°

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The Beauty Tech Group Ltd.