



neonatal
solutions

Directions for Use

Enteral Accessories

ENFit®

R_x Only

Not made with natural rubber latex.

Not made with Diethylhexylphthalate (DEHP).

WARNING:

- The ENFit Enteral Accessories are not intended for intravascular or intravenous use.
- This device has the potential to misconnect with small bore connectors of other healthcare applications.

CAUTIONS:

- Contents are sterile in an unopened, undamaged package. Do not use this device if opened or damaged.
- Do not re-sterilize.

INDICATIONS FOR USE:

- The ENFit Enteral Accessories are intended to be attached to an enteral extension set, nasogastric feeding tube and/or enteral syringe. The ENFit Enteral Accessories are disposable and single use only.

RECOMMENDED PROCEDURE:

- Use clean technique.
- Ensure this device is only connected to an enteral port and not to an I.V. set.

ENTERAL SET REMOVAL:

- Remove the enteral accessory and discard per hospital protocol

 Medela LLC
1101 Corporate Drive,
McHenry IL 60050

Medela wordmark and logo is registered in the U.S. Patent and Trademark Office. ENFit® is a federally registered trademark of GEDSA and is used with permission.

1548594 A 0617 © 2017 Medela DFU-045