

vive
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

WHEELCHAIR BAG STANDARD

LVA3098



To see all of the FAQs in one place,
visit vhealth.link/f4d156

vive
vivehealth.com

QUICK START GUIDE

290525

INTENDED USE

The Vive Wheelchair Bag provides convenient storage for personal items on wheelchairs and transport chairs.

APPLICATION INSTRUCTIONS

1. ATTACHING THE BAG

- Slip the two straps over the wheelchair's push handles.
- Position the bag so the flap faces outward.



2. USING THE BAG

- Lift the flap to place items in the storage compartment.
- Close the flap and secure with the hook and loop fastener.



CARE INSTRUCTIONS

- Machine wash in cold water
- Air dry only
- Do not bleach or iron

WARNINGS ⚠

- Maximum weight capacity: 25 lbs
- Not waterproof - keep away from rain and moisture.
- Do not use to carry medical devices that require secure positioning.
- This product is not intended to prevent, cure, or treat any medical condition.

WARRANTY

Your Vive Wheelchair Bag has a 1-year warranty from date of purchase. This warranty covers manufacturing defects and workmanship. Normal wear and tear is not covered.